

GLOBE Electronic Precision Balance



CS3001 3000g x 0.1g

CS302 300g x 0.01g

CS602 600g x 0.01g

Standard Features:

- AC/DC can be used
- 11 selectable weighing units
- CS302 & CS602 have 2 additional fabric weighing units: oz/yd² & g/m²
- Protection screw designed for safety transportation
- CE Approved
- RS-232 Uni-direction for data output
- Wind-Shield included (Removable)
- Weighing, counting and percentage modes
- Zero tracking function
- Two overload protection systems
- LCD display backlight
- Temperature compensation
- 15 range of stability
- Full capacity taring

Specification:

Model	CS302	CS602	CS3001
Capacity	300g	600g	3000g
Readability	0.01g	0.01g	0.1g
Weighting Units ***	g, ct, lb, oz, oz/yd ² *, g/m ² *, ozt, dwt, MM(Momme), 3 Taels, t(Tola)		
Platform Size	Φ 116 mm/ Φ 4.5"		124 x 144 mm /4.8"x 5.6"
Power	110/220V AC Adaptor** (Supplied) OR 4 x UM-2 Battery (Not Supplied)		
Dimension	226(L) x 175(W) x 136(H) mm / 8.8" (L) x 6.8" (W) x 5.3" (H)		
Weight	Approx. 1.1kg/2.42lb (Included Wind-Shield)		
Working Temperature	-10 ~ 40°C		

* For Model CS302, CS602 only

** Please specify voltage when order

***The balance permits you to specify which weighing units and operating modes will be enabled for use during operation

GLOBE Economy Electronic Balance



CS202 200g x 0.01g
EB302 300g x 0.01g

CS2001 2000g x 0.1g

CS5000 5000g x 1g

Standard Features:

- Low Price
- AC/DC can be used
- Auto Zero Tracking
- Net weight/Stability Indication
- Auto shut off
- Model CS202 & EB302 included Wind shield
- Low Battery indication
- Auto backlight
- Large LCD (6 digits 15mm high)

Specification:

Model	CS202	EB302	CS2001	CS5000
Capacity	200g	300g	2000g	5000g
Readability	0.01g	0.01g	0.1g	1g
Weighting Units ***	g, g/m2*, oz, oz/yd2*		g, ct, oz, dwt	
Platform Size	Φ 100 mm/ Φ 4"		145 x 145 mm / 5.6" x 5.6"	
Wind Shield	Supplied		Not Supplied	
Power	110/220V AC Adaptor** (Supplied) OR 6 x AA Battery (Not Supplied)			
Dimension	Main Body Size : 145 x 210 x 40mm / 5.6" x 8.2" x 1.6" Wind Shield Size : 145 x 145 x 120mm / 5.6" x 5.6" x 4.6" (CS202 & EB302 ONLY)			
Weight	Approx. 1kg/2.2lb			
Working Temperature	0 ~ 40°C			

* Please specify voltage when order