

Blood Pressure Set, Model PA-130, this standard aneroid sphygmomanometer is a precise blood pressure instrument with a 300 mm aneroid gauge. Complete with a velcro-type cuff inflation system and rubber bulb with release valve. Supplied in a zippered carrying case.



PA-130



PA-130A

Mercury Sphygmomanometer, Model PA-130A, this is a standrd mercury column blood pressure instrument with a 300 mm scale. The glass tube is made of heavy glass to reduce chances of breakage. The instrument features a cuff, a special shutoff valve which prevents the loss of mercury, and a sturdy metal case.


Model PA-SS-1
Single head

PA-SS-2
Double heads



PA-SS-1 & PA-SS-2

Stethoscope, this sensitive stethoscope has a 4cm chestpiece with a raised stem, 56cm PVC Y-tube, and adjustable chrome-plated brass binaurals, white ear tips.



PA-130B

Aneroid Sphygmomanometer with built-in stethoscope, Model PA-130B, just clip the cuff over your arm and secure fastener with one hand. Listen with the built-in stethoscope and take the reading to 300mm with an aneroid manometer. Set comes with vinyl zipper case and includes instructions and charts.



PA-DS-2

Teaching Stethoscope, Model PA-DS-2, designed so both you and a student can monitor a sound simultaneously, this instrument has a dual head with two sets of chrome-plated brass binaurals and two 56 cm Y-tubes. A metal piece connects the Y-tubes to the 25 cm extension tube and chestpiece.



PA-FS-1

Stethoscope Ford Type, Model PA-FS-1, this ford type stethoscope has a deep 2.5 cm bell, dual stem chestpiece, and adjustable chrome-plated brass binaurals, sensitive acoustic response, making this an excellent instrument for listening to heart and chest sounds.