


Two pans student balance, **Model PA-250**, ideal for student use with weight set, capacity : 200gm., sensitivity : 20mg.,



**Three beams balance**, single pan, flat stainless steel plate pan 150mm diameter, capacity 2610g., sensitivity 0.1g., with three beams graduated 0 to 500g x 100g., 0 to 100g x 10g and 0 to 10g x 0.1g. Masses are unnecessary for weighing up to 610 g. Complete with additional supplementary masses (set of 3 pcs - 2 x 1000g and 1 x 500g) to hang on the end of the beam increasing the capacity up to 2610g. Zero adjustment, magnetically damped beam movement ensures rapid weighing. (balance and weight sets can be sold separately).

**Model PA-251** without additional ungraduated tare beam.

**Model PA-251A** with additional ungraduated tare beam, capacity 225g.,

**Balance, single pan, Model PA-252**, capacity 311g., sensitivity 0.01g., with four beams, each provided with a non-removable slider and graduated 0 to 200g x 100g, 0 to 100g x 10g, 0 to 10g x 1g and 0 to 1 g x 10mg respectively. Additional masses are not required. With zero adjustment and stainless steel pan 90mm diameter with two pouring spouts and handle. Height of stirrup 175mm. Magnetically damped beam movement ensures rapid weighing. A platform housed in the base can be used for density determinations.



PA-252



BR-9011

**Bathroom Scale, Model BR-9011**, the easy-to-read dial on this bathroom scale is calibrated in kilograms and pounds, and registers to 125 kg (275 pounds). The ruggedly constructed scale has a heavy metal housing and black rubber foot mat.