

Laboratory Porcelain Wares

The Chinese Laboratory Porcelain Ware are made of high purity raw material. The Products thus produced are uniform in quality, compact in structure, impervious to gas infiltration, good in polish and resistant to acids & alkalis, with sufficient mechanical strength, the porcelain can withstand sudden temperature changes without cracks, explosion or deformation. Under normal condition, the glazed pieces can sustain a temperature of 1050 deg C while the unglazed ones a temperature of 1350 deg C. All products are glazed inside and outside except outside bottom.



PW-1 to PW-7



PW-8 to PW-10



PW-11 to PW-18



PW-24 to PW-24A



PW-28 to PW-32

Model

- PW-1** Porcelain crucible. high form. with cover (lid)., cap. 10 ml. rim dia. 30 mm. Ht: 28 mm.
- PW-2** - ditto - but cap. 15 ml. rim dia. 35 mm. Ht: 30 mm.
- PW-3** - ditto - but cap. 18 ml. rim dia. 37 mm. Ht: 31 mm.
- PW-4** - ditto - but cap. 25 ml. rim dia. 40 mm. Ht: 36 mm.
- PW-5** - ditto - but cap. 30 ml. rim dia. 42 mm. Ht: 41 mm.
- PW-6** - ditto - but cap. 50 ml. rim dia. 53 mm. Ht: 46 mm.
- PW-7** - ditto - but cap. 100 ml. rim dia. 63 mm. Ht: 59 mm.
- PW-8** Porcelain crucible. low wide form. with cover (lid)., cap. 15 ml. rim dia. 43 mm. Ht: 23 mm.
- PW-9** - ditto - cap. 25 ml. rim dia. 45 mm. Ht: 25 mm.
- PW-10** - ditto - cap. 30 ml. rim dia. 51 mm. Ht: 30 mm.
- PW-11** Porcelain Basin (Evaporating Dish). round bottom, withspout. cap. 35 ml. rim diam. 60 mm. Ht: 25 mm.
- PW-12** - ditto - cap. 50 ml. rim dia. 70 mm. Ht: 28 mm.
- PW-13** - ditto - cap. 75 ml. rim dia. 82 mm. Ht: 32 mm.
- PW-14** - ditto - cap. 100 ml. rim dia. 90 mm. Ht: 35 mm.
- PW-15** - ditto - cap. 125 ml. rim dia. 98 mm. Ht: 38 mm.
- PW-16** - ditto - cap. 150 ml. rim dia. 104 mm. Ht: 40 mm.
- PW-17** - ditto - cap. 200 ml. rim dia. 110 mm. Ht: 44 mm.
- PW-17A** - ditto - cap. 250 ml. rim dia. 118 mm. Ht: 48 mm.
- PW-18** - ditto - cap. 270 ml. rim dia. 120 mm. Ht: 49 mm.
- PW-24** Gooch Crucible with cover (lid) & perforated disc., cap. 25 ml.
- PW-24A** - ditto - cap. 30 ml.
- PW-28** Porcelain Buchner Funnel. with perforated plate, 50 ml. int. dia. 51 mm.
- PW-29** - ditto - but 120 ml. int. dia. 70 mm.
- PW-30** - ditto - but 200 ml. int. dia. 90 mm.
- PW-31** - ditto - but 400 ml. int. dia. 110 mm.
- PW-31A** - ditto - but 600 ml, int. dia., 125mm.
- PW-32** - ditto - but 1000 ml. int dia. 150 mm.



PW-25, PW-33A to PW-35



PW-38 & PW-38A



PW-36 to PW-37



PW-45 to PW47



PW-39 to PW-41



PW-48

New

PW-49



New

Model

- PW-33A** Porcelain Mortar & Pestle. with spout. 60 ml. rim dia. 80 mm Ht: 40 mm.
- PW-33** - ditto - 100 ml. rim dia. 100 mm. Ht: 45 mm.
- PW-34** - ditto - 300 ml. rim dia. 130 mm. Ht: 62 mm.
- PW-35** - ditto - 500 ml. rim dia. 160 mm. Ht: 74 mm.
- PW-25** - ditto - 1000ml. rim dia. 216mm, Ht : 90mm
- PW-36** Porcelain combustion boat. Unglazed, length 72 mm. Ht: 9 mm. rim width: 11 mm.
- PW-37** length 88 mm. Ht: 10 mm. rim width: 13 mm.
- PW-38** Porcelain Spotting Plate. with 12 depressions, size: 115 x 90 x 12 mm. (L x W x Thickness).
- PW-38A** - ditto -, with 6 depression, size : 85 x 54 x 11 mm,
- PW-39** Pipe Clay Triangle, with iron wire, length of side : 50 mm.
- PW-40** - ditto -, length of side : 63 mm.
- PW-41** - ditto -, length of side : 75 mm.
- PW-45** Desiccator porcelain plate, with numerous 10mm (approx.) perforations, large hole in center, glazed on one side, dia. 140mm.
- PW-46** Desiccator porcelain plate, with numerous 10mm (approx.) perforations, large hole in center, glazed on one side, dia. 190mm.
- PW-47** Desiccator porcelain plate, with numerous 10mm (approx.) perforations, large hole in center, glazed on one side, dia. 230mm
- PW-48** Beehive Shelves, diam. 100mm
- PW-49** Acid tray, diam. 100mm, H: 20mm..