



PA-250

Two pans student balance, **Model PA-250**, ideal for student use with weight set, capacity : 200gm., sensitivity : 20mg.,




PA-252

Balance, single pan, Model PA-252, capacity 311g., sensitivity 0.01g., with four beams, each provided with a non-removable slider and graduated 0 to 200g x 100g, 0 to 100g x 10g, 0 to 10g x 1g and 0 to 1 g x 10mg respectively. Additional masses are not required. With zero adjustment and stainless steel pan 90mm diameter with two pouring spouts and handle. Height of stirrup 175mm. Magnetically damped beam movement ensures rapid weighing. A platform housed in the base can be used for density determinations.



PA-322 and PA-323

“BONSO” Electronic Balance, a range of low cost high quality top pan balances. They have a large easy to read LCD display and big stainless steel weighing pan. Operation is by two buttons.



PA-251 & PA-251A

Three beams balance, single pan, flat stainless steel plate pan 150mm diameter, capacity 2610g., sensitivity 0.1g., with three beams graduated 0 to 500g x 100g., 0 to 100g x 10g and 0 to 10g x 0.1g. Masses are unnecessary for weighing up to 610 g. Complete with additional supplementary masses (set of 3 pcs - 2 x 1000g and 1 x 500g) to hang on the end of the beam increasing the capacity up to 2610g. Zero adjustment, magnetically damped beam movement ensures rapid weighing. (balance and weight sets can be sold separately).

Model PA-251 without additional ungraduated tare beam, capacity 225g.,

Model PA-251A with additional ungraduated tare beam, capacity 225g.,



BR-9011

Bathroom Scale, Model BR-9011, the easy-to-read dial on this bathroom scale is calibrated in kilograms and pounds, and registers to 125 kg (275 pounds). The ruggedly constructed scale has a heavy metal housing and black rubber foot mat.

Model	PA-322	PA-323
Capacity	2000g	5000g
Graduation	1g	1g
Weight modes	g/lb-oz	gram only
Power supply (optional)	1 x 9 volt alkaline battery or UL/CSA approved AC adaptor (both are not included)	

Display size	11 mm x 4 digits
Pan size :	150 x 170 mm
Dimensions :	230 x 165 x 80 mm
features	Auto power off, tare function, includes stainless steel tray top