

## CREST POWERSONIC™ BENCHTOP CLEANERS SPECIFICATIONS

Model	Inside Tank Dimensions L x W x H inches centimeters	Overall Dimensions L x W x H inches centimeters	Tank Capacity gallon liter	Ship Wt. lbs kilos	Line Input @117V amps (rated)			Sonic Power Avg. Wattage
					T	HT	D	
CP 200	6 x 5.25 x 4 148 x 134 x 100	7 x 6.5 x 10.5 175 x 164 x 270	0.5 1.75	6 2.7	0.8	1.5	N/A	55
CP 230	9.5 x 5.25 x 4 237 x 134 x 100	11.5 x 6.375 x 9.25 265 x 162 x 295	0.75 2.8	9 4.1	1.2	2.6	2.6	80
CP 360	9.5 x 5.25 x 6 237 x 134 x 150	11.5 x 6.375 x 11.625 265 x 162 x 295	1 3.8	10.8 4.9	1.5	2.7	2.7	100
CP 500	11.75 x 6 x 6 297 x 148 x 150	12.75 x 7 x 11.625 325 x 176 x 295	1.5 5.4	12.8 5.8	1.7	3	3	120
CP 1200	19.5 x 5.25 x 6 495 x 134 x 150	20.25 x 6 x 10.625 515 x 150 x 270	2.5 9.5	21 3.5	2.6	7.8	7.8	200
CP 1100	11.75 x 9.5 x 8 297 x 237 x 200	12.75 x 10.5 x 13 325 x 265 x 335	3.25 12.3	22 10	2.6	7.8	7.8	200
CP 2600	19.5 x 11.75 x 8 495 x 297 x 200	21 x 12.75 x 14.5 530 x 325 x 365	7 26.1	38 17.3	3.8	12	12	300

**NOTES:**

- All models are available in three (3) configurations:
  - T** 0-30 minute timer
  - HT** 0-30 minute times and thermostatically controlled heat ambient to 80°C (175°F)
  - \*D** Digital control of heat, time (0-99 minutes), power and degas functions
- Frequency 42kHz – 45kHz
- All models meet FCC requirements
- All models are ETL/CSA approved
- All models are available in 117 volts and 230 volts

\* Digital controls not available with CP200 model



CP 1200HT with SSMB 1200-DH Mesh Basket with Drain Positioning Handles



**香港科學儀器社有限公司**  
Che Scientific Co. (Hong Kong) Ltd.

**Hong Kong Head office:**

Unit D, 18/F, Roxy Industrial Centre, 58-66 Tai Lin Pai Road, Kwai Chung, Hong Kong

TEL: (852) 2481 1323 FAX: (852) 2418 1302

E-mail: [sales@chescientific.com](mailto:sales@chescientific.com) Website: [www.chescientific.com](http://www.chescientific.com)



**浙江科通儀器有限公司**  
Zhejiang Fortune Scientific Co. Ltd.

TEL: (86 571) 8729 7500 FAX: (86 571) 8729 7511

E-mail: [sales@fortunescientific.com.cn](mailto:sales@fortunescientific.com.cn)

Website: [www.fortunescientific.com.cn](http://www.fortunescientific.com.cn)

**Guangzhou Branch:**

TEL: (86 20) 8130 7388 FAX: (86 20) 8108 4925

E-mail: [guangzhou@fortunescientific.com.cn](mailto:guangzhou@fortunescientific.com.cn)