


WEIFO® "WEIFO" Electronic Auto-Drying Box Made in Taiwan

WEIFO dry box condenses moisture in the box and exhaust the moisture to outside with the T.E.COOLING system. Full auto humidity control provides the best storage stage. Perfectly solve the problems of damp, mould, rust, metamorphic and aging.

Practical Standard Series Electronic Auto-Drying Box



Economy Series Electronic Auto-Drying Box



Model	Practical Standard Series (Digital humidity indicator, 20%~60% RH setting)				Economy Series (Bar type humidity indicator, 30/40/50/60% RH setting)	
	DRY-60	DRY-70	DRY-100	DRY-205	KD-70	KD-112
External Dim. (WxDxHmm)	400x410x560	400x410x660	440x410x890	400x530x1280	400x410x560	400x410x890
Internal Dim. (WxDxHmm)	398x355x480	398x355x580	398x355x820	398x475x1210	398x355x480	398x355x820
Capacity	68L	82L	116L	228L	70L	116L
Door	SINGLE	SINGLE	DUAL	DUAL	SINGLE	SINGLE
Shelf	2 PCS	2 PCS	3 PCS	5 PCS	2 PCS	3 PCS
Power Cons. (Max/Average)	12W/4W	12W/5W	12W/6W	22W/10W	12W/4W	12W/5W
Weight	15kg	17kg	20kg	30kg	15kg	20kg

Remark: No temperature indicator for practical standard series and economy series.
SuperDry low humidity type is available. (Humidity range: 10%~60% RH)

Professional Series Electronic Auto-Drying Box



PROFESSIONAL SERIES (With humidity and temperature display)							
Model	DRY-125	DRY-300	DRY-420	DRY-520	DRY-610	DRY-530	DRY-1360
External Dim. (WxDxHmm)	440x440x930	890x480x930	1180x480x1020	1180x570x1020	1400x600x100 With castor	920x660x1190 With castor	1200x730x1900 With castor
Internal Dim. (WxDxHmm)	438x385x850	888x425x850	1178x425x900	1178x515x900	1398x545x850	888x610x1070	1198x675x1740
Capacity	141L	320L	450L	545L	646L	578L	1402L
Door	SINGLE	DUAL	DUAL	DUAL	DUAL	4 PCS Doors (Front & Back)	6 PCS Doors
Shelf	3 PCS	3 PCS	3 PCS	3 PCS	3 PCS	6 PCS	5 PCS
Power Cons. (Max/Average)	12W/6W	22W/10W	32W/12W	32W/14W	32W/16W	32W/16W	52W/25W
Weight	25kg	43kg	55kg	60kg	68kg	70kg	120kg

Remark: SuperDry low humidity type is available. (Humidity range: 10%~60% RH)



香港科學儀器社有限公司
Che Scientific Co. (Hong Kong) Ltd.

Hong Kong Head office:

Unit D, 18/F, Roxy Industrial Centre, 58-66 Tai Lin Pai Road, Kwai Chung, Hong Kong

TEL: (852) 2481 1323 FAX: (852) 2418 1302

E-mail: sales@chescientific.com Website: www.chescientific.com



浙江科通儀器有限公司
Zhejiang Fortune Scientific Co. Ltd.

TEL: (86 571) 8729 7500 FAX: (86 571) 8729 7511

E-mail: sales@fortunescientific.com.cn

Website: www.fortunescientific.com.cn

Guangzhou Branch:

TEL: (86 20) 8130 7388 FAX: (86 20) 8108 4925

E-mail: guangzhou@fortunescientific.com.cn