



“格拉斯科”牌 科学仪器  
“Glas-Col” Science Equipment

美国  
USA

## TM Heating Mantles 电热罩

- Rigid aluminum housing provides strength and durability while supporting the weight of the vessel
- Fabric interior to softly nest glass vessels and reduce the chance of thermal shock
- Adaptable to most vessels including larger sizes
- Withstands 450°C internal operating temperature (650°C for Series STM), enough power for most applications
- Optional temperature controller available (Sold separately)



控制器需另购

Specification 规格							
Catalog Number 230 volt	Flask Capacity (ml)	Max. Flask Diameter (mm)	Watts	Inside Depth (in)	Outside Diameter (in)	Outside Height (in)	Product Weight (lbs)
100B TM97	100	60	80	1.19	6.25	4.13	2.2
100B TM99	125	70	80	1.38	6.25	4.13	2.2
100B TM101	200	76	100	1.5	6.25	4.13	2.2
100B TM103	250	83	180	1.63	6.25	4.75	2.5
100B TM105	300	86	180	1.69	6.25	4.75	2.5
100B TM107	500	101	270	2	6.25	5	2.5
100B TM109	1,000	137	380	2.56	7.5	5	3
100B TM111	2,000	170	500	3.38	10	6	4.7
100B TM113	3,000	183	500	3.6	10	6	4.7
100B TM115	5,000	220	600	4.31	12	7	6.2
100B TM117	12,000	293	2@650	5.77	16.5	9	15
100B TM119	22,000	347	2@770	6.81	19	10	18.8
100B TM121	50,000	456	3@1,000	8.98	24	13	41
100B TM122	72,000	522	2@2,000	10.19	26	14	44

## COMBO Mantles 电热罩

- Six-place for 100-300 ml flasks
- Lower profile for space-saving convenience
- Revolutionary heating element container system makes replacing burned-out elements a simple, five minute task
- Element containers act as spill containment chambers to capture chemicals from spillage or boil-over
- Clear anodized aluminum cabinet with black PTFE resin coated top for chemical resistance
- Available in 450°C and 650°C version for greater versatility
- Optional temperature controller available (Sold separately)



Specification 规格							
Catalog Number 230 volt	Flask Capacity (ml)	Max. Flask Diameter (in)	Max. temperature (°C)	Bottom	Rating per position (W)	Replacement Element No.	Product Weight (lbs)
100D RX30424	250/300	3.41	450	Round	125	100D RX30424E	19.0
100D RJ30424	250/300	3.41	650	Round	210	100D RJ30424E	19.0

## 3D Shaker 三维摇动器

- Mechanically shakes glassware for consistent, repeatable results.
- Handles total loads up to 40 lbs without needing to balance load.
- "Lazy Susan" design allows 360° rotating of platform for easy loading and unloading.
- Platform locks in place for added safety.
- Accommodates up to four 2-liter separatory funnels.
- Glas-Col's automatic venting separatory funnels available for added safety



Specification 规格			
Catalog Number 230 volt	Description	RPM	Weight (lbs)
099A VS21024	Shaker Base with four 1-liter separatory funnel holders	10-170	165
099A VS22024	Shaker Base with four 2-liter separatory funnel holders	10-170	165

香港科學儀器社有限公司  
Che Scientific Co. (Hong Kong) Ltd.  
Tel: (852) 2481 1323  
Fax: (852) 2418 1302

浙江科通儀器有限公司(杭州)  
Zhejiang Fortune Scientific Co. Ltd. (Hangzhou)  
Tel: (86 571) 8729 7500  
Fax: (86 571) 8729 7511

浙江科通儀器有限公司(廣州)  
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)  
Tel: (86 20) 8108 9804  
Fax: (86 20) 8108 4925

上海香科儀器貿易有限公司  
Shanghai Xiang Ke Scientific Trading Co. Ltd.  
Tel: (86 21) 6248 6170  
Fax: (86 21) 6248 0647