

PYREX®

BOECO
Germany



实验室玻璃设备
Laboratory Glassware

KIMBLE
KONTES

JAYTEC

Low / Tall Form Beakers 矮身/高身烧杯

Capacity 容量	Low Form 矮身	Tall Form 高身
10 mL	✓	
25 mL	✓	
50 mL	✓	✓
100 mL	✓	✓
150 mL	✓	✓
250 mL	✓	✓
400 mL	✓	✓
600 mL	✓	✓
800 mL	✓	✓
1 L	✓	✓
2 L	✓	✓
3 L	✓	✓
5 L	✓	✓



Borosilicate Glass

Wide Mouth Reagent Bottles 阔口试剂瓶

Capacity 容量
25 mL
50 mL
100 mL
250 mL
500 mL
1 L
2 L
5 L
10 L
20 L



Borosilicate Glass

Blue screw cap autoclavable at 121°C
蓝色螺丝胶盖可于 121°C 高温消毒

Conical Flasks 三角烧瓶

Capacity 容量	N/M 窄口	W/M 阔口
25 mL	✓	
50 mL	✓	✓
100 mL	✓	✓
200 mL	✓	✓
250 mL	✓	✓
500 mL	✓	✓
1 L	✓	✓
2 L	✓	✓
3 L	✓	✓
5 L	✓	✓



Borosilicate Glass

Round / Flat Bottom Flasks 圆底/平底烧瓶

Capacity 容量
50 mL
100 mL
250 mL
500 mL
1 L
2 L
4 L
6 L
10 L
20 L



Borosilicate Glass

PYREX® Filtering Flask 过滤瓶

Heavy wall use for vacuum filtration
特厚 PYREX® 玻璃适当真空过滤使用


Capacity 容量	Neck ID 瓶口内径
100 mL	24 mm
250 mL	35 mm
500 mL	35 mm
1 L	45 mm
2 L	60 mm
5 L	68 mm



Desiccators 干燥器

Vacuum model with stopcock
真空干燥器附有活塞

Inner dia. 内径
150 mm
200 mm
250 mm
300 mm




Standard

Vacuum

Measuring Cylinder 量筒

Capacity 容量	Graduations 最少刻度	Tolerance 误差 (±mL)	
		Class A	Class B
5 mL	0.1 mL	0.05	0.1
10 mL	0.2 mL	0.10	0.2
25 mL	0.5 mL	0.25	0.5
50 mL	1.0 mL	0.50	1.0
100 mL	1.0 mL	0.50	1.0
250 mL	2.0 mL	1.0	2.0
500 mL	5.0 mL	2.5	5.0
1 L	10.0 mL	5.0	10.0
2 L	20.0 mL	10.0	20.0



Volumetric Flask 定量瓶

Capacity 容量	Tolerance 误差 (±mL)	
	Class A	Class B
5 mL	0.025	0.05
10 mL	0.025	0.05
25 mL	0.04	0.08
50 mL	0.06	0.12
100 mL	0.10	0.20
250 mL	0.15	0.30
500 mL	0.25	0.50
1 L	0.40	0.80
2 L	0.60	1.20

