



## “奥斯摩克斯”牌 过滤膜及过滤装置 “OSMONICS” Filter Devices, Membranes

美国制造  
Made in USA

### Gravi-Seal Analytical Funnel 分析过滤漏斗

- 只需单手操作，稳定易用
- 47mm 直径过滤膜使用
- 容量 350ml，每 50ml 有刻度
- 无需夹架和磁铁等消耗品
- 宽阔的过滤膜提取点，方便快速取出过滤膜



物料	完整套装 (1PC)	特惠套装 (3PC)
Polysulfone	#PFC0004701	#PFC0004701
Polypropylene	#PPC0004701	#PFC0004701

### MicronSep 1+PAC White Filter Membrane 白过滤膜 (独立包装)

- Designed for the recovery and retention of E.coli bacteria in water / wastewater analysis and other microbiological tests.
- The filters are certified to meet specifications list in APHA Standard Methods.
- Each 1+PAC includes a 47mm presterilized MicronSep nitrocellulose filter.



滤孔大小 Pore Size	表面 Surface	吸收剂垫 Absorbent Pad	产品编号 Catalog No.	
0.22 μm	Grid	Yes	#E02WG047S2 100/PK	#E02WG047S5 1000/PK
0.22 μm	Grid	No	#E02WG047S1 200/PK	#E02WG047S4 1000/PK
0.45 μm	Grid	Yes	#E04WG047S2 100/PK	#E04WG047S5 1000/PK
0.45 μm	Grid	No	#E04WG047S1 200/PK	#E04WG047S4 1000/PK
0.70 μm	Grid	Yes	#E07WG047S2 100/PK	#E07WG047S5 1000/PK
0.70 μm	Grid	No	#E07WG047S1 200/PK	#E07WG047S4 1000/PK


### MicronSep 1+PAC Black Filter Membrane 黑过滤膜 (独立包装)

- Designed for quality assurance testing in food and beverage analysis.
- Provide a superb contrast for early and accurate colony counting for total bacteria, yeast and mold organisms.
- The filters are certified to meet specifications list in APHA Standard Methods.
- Each 1+PAC includes a 47mm presterilized MicronSep black nitrocellulose filter.



滤孔大小 Pore Size	表面 Surface	吸收剂垫 Absorbent Pad	产品编号 Catalog No.	
0.22 μm	Grid	Yes	#E02WG047S2 100/PK	#E02WG047S5 1000/PK
0.22 μm	Grid	No	#E02WG047S1 200/PK	#E02WG047S4 1000/PK
0.45 μm	Grid	Yes	#E04WG047S2 100/PK	#E04WG047S5 1000/PK
0.45 μm	Grid	No	#E04WG047S1 200/PK	#E04WG047S4 1000/PK
0.70 μm	Grid	Yes	#E07WG047S2 100/PK	#E07WG047S5 1000/PK
0.70 μm	Grid	No	#E07WG047S1 200/PK	#E07WG047S4 1000/PK

### AcetatePlus Cellulose Acetate Membrane Filters 有载体醋酸纤维素过滤膜



- Protein and enzyme filtration, sterilization
- Biological fluid filtration sterilization
- Tissue culture media sterilization

直径 Diameter	滤孔大小 Pore Size	表面 Surface	产品编号 Catalog No.
47mm	0.22 μm	Plain	#A02SP04700 100/PK
47mm	0.45 μm	Plain	#A04SP04700 100/PK

### PreSep G40 Glass Fiber Prefilters 玻璃纤维前过滤膜 (No binder)



- For solvent filtration, gravimetric analysis, contamination analysis of wastewater, etc.
- Suitable for higher temperatures
- Autoclavable

直径 Diameter	滤孔大小 Pore Size	产品编号 Catalog No.
47mm	1.0 μm	#G40WP04705 500/PK
90mm	1.0 μm	#G40WP09000 100/PK




## “奥斯摩克斯”牌 过滤膜及过滤装置 “OSMONICS” Filter Devices, Membranes

美国制造  
Made in USA

### Magna Nylon Membranes

- Sterilization, clarification of aqueous and organic solvent solutions
- HPLC sample preparation and solvent filtration

直径 Diameter	滤孔大小 Pore Size	表面 Surface	产品编号 Catalog No.
47mm	0.22 μm	Plain	#R02SP04700 100/PK
47mm	0.45 μm	Plain	#R04SP04700 100/PK



### Magna PTFE (Teflon) Laminated Membranes

- Filtration of strong acids and aggressive solutions
- Venting applications

直径 Diameter	滤孔大小 Pore Size	厚度 Thickness	产品编号 Catalog No.
47mm	0.22 μm	175 μm	#F02LP04700 100/PK
47mm	0.45 μm	175 μm	#F05LP04700 100/PK
47mm	1.00 μm	175 μm	#F10LP04700 100/PK



### Poretics® Polycarbonate (PCTE) Membrane Filters

- Water analysis, oceanography
- Surface capture for analysis
- Low protein binding applications

直径 Diameter	滤孔大小 Pore Size	产品编号 Catalog No.
47mm	1.0 μm	#K10CP04700 100/PK
47mm	2.0 μm	#K20CP04700 100/PK
47mm	3.0 μm	#K30CP04700 100/PK
47mm	5.0 μm	#K50CP04700 100/PK
47mm	8.0 μm	#K80CP04700 100/PK



### Cameo Acetate Syringe Filter, Low Protein Binding

- Sterilize biological fluids, serum or media additives
- Sample preparation of aqueous solutions
- HPLC sample preparation

直径 Diameter	滤孔大小 Pore Size	消毒 Sterile	产品编号 Catalog No.
25mm	0.22 μm	Yes	#DDA02025S0 50/PK
25mm	0.45 μm	Yes	#DDA04025S0 50/PK
30mm	0.22 μm	No	#DDA02030S0 50/PK
30mm	0.45 μm	No	#DDA04030S0 50/PK



### Cameo Nylon Syringe Filter

- Sterilize aqueous or/and organic samples
- Sterilize and clarification biological fluids, buffer solution and tissue culture media
- HPLC aqueous and organic sample filtration

直径 Diameter	滤孔大小 Pore Size	消毒 Sterile	产品编号 Catalog No.
25mm	0.22 μm	Yes	#DDR02025S0 50/PK
25mm	0.45 μm	Yes	#DDR04025S0 50/PK
30mm	0.22 μm	No	#DDR02T30S0 50/PK
30mm	0.45 μm	No	#DDR04T30S0 50/PK



### Cameo Teflon Syringe Filter

- Aggressive solvent, acids and bases, organic-based HPLC samples
- Venting applications, sterilize air, gas, aerosols

直径 Diameter	滤孔大小 Pore Size	消毒 Sterile	产品编号 Catalog No.
25mm	0.22 μm	No	#DDF02T1750 50/PK
25mm	0.45 μm	No	#DDF04T1750 50/PK
30mm	0.22 μm	No	#DDF02T3050 50/PK
30mm	0.45 μm	No	#DDF04T3050 50/PK

