

Balance and Weighing Scales 天平及磅称

TROEMNER

“特朗纳”牌 标准砝码 “TROEMNER” Standard Weights

美国制造
Made in USA

OIML Standard Weights 标准砝码 (OIML – Organization Internationale de Metrologie Legale)



Tolerances 公差 (in \pm mg)				
Calibration Tolerances International Recommendation OIML R111				
Nominal Value	E1 Class	E2 Class	F1 Class	F2 Class
1mg	0.002	0.006	0.020	0.06
2mg	0.002	0.006	0.020	0.06
5mg	0.002	0.006	0.020	0.06
10mg	0.002	0.008	0.025	0.08
20mg	0.003	0.010	0.03	0.10
50mg	0.004	0.012	0.04	0.12
100mg	0.005	0.015	0.05	0.15
200mg	0.006	0.020	0.06	0.20
500mg	0.008	0.025	0.08	0.25
1g	0.010	0.030	0.10	0.3
2g	0.012	0.040	0.12	0.4
5g	0.015	0.050	0.15	0.5
10g	0.020	0.060	0.20	0.6
20g	0.025	0.080	0.25	0.8
50g	0.030	0.10	0.30	1.0
100g	0.05	0.15	0.5	1.5
200g	0.10	0.30	1.0	3.0
500g	0.25	0.75	2.5	7.5
1kg	0.5	1.5	5	15
2kg	1.0	3.0	10	30
5kg	2.5	7.5	25	75
10kg	5	15	50	150
20kg	10	30	100	300
50kg	25	75	250	750

Order Information 订购资料

Individual Standard Weights 单一标准砝码					Standard Weight Sets 标准砝码套装				
Nominal Value	E1 Class	E2 Class	F1 Class	F2 Class	Nominal Value	E1 Class	E2 Class	F1 Class	F2 Class
1mg	7537-E1P	7537-E2P	7537-F1P	7537-F2P	500mg – 1mg	7440-E1P	7440-E2P	7440-F1P	7440-F2P
2mg	7536-E1P	7536-E2P	7536-F1P	7536-F2P	5g – 1mg	7441-E1P	7441-E2P	7441-F1P	7441-F2P
5mg	7534-E1P	7534-E2P	7534-F1P	7534-F2P	50g – 1mg	7433-E1P	7433-E2P	7433-F1P	7433-F2P
10mg	7533-E1P	7533-E2P	7533-F1P	7533-F2P	50g – 1g	7429-E1P	7429-E2P	7429-F1P	7429-F2P
20mg	7532-E1P	7532-E2P	7532-F1P	7532-F2P	100g – 1mg	7428-E1P	7428-E2P	7428-F1P	7428-F2P
50mg	7530-E1P	7530-E2P	7530-F1P	7530-F2P	100g – 1g	7424-E1P	7424-E2P	7424-F1P	7424-F2P
100mg	7529-E1P	7529-E2P	7529-F1P	7529-F2P	200g – 1mg	7423-E1P	7423-E2P	7423-F1P	7423-F2P
200mg	7528-E1P	7528-E2P	7528-F1P	7528-F2P	200g – 1g	7422-E1P	7422-E2P	7422-F1P	7422-F2P
500mg	7526-E1P	7526-E2P	7526-F1P	7526-F2P	500g – 1mg	7419-E1	7419-E2	7419-F1	7419-F2
1g	7525-E1P	7525-E2P	7525-F1P	7525-F2P	500g – 1g	7418-E1	7418-E2	7418-F1	7418-F2
2g	7524-E1P	7524-E2P	7524-F1P	7524-F2P	1kg – 1mg	7417-E1	7417-E2	7417-F1	7417-F2
5g	7522-E1P	7522-E2P	7522-F1P	7522-F2P	1kg – 1g	7416-E1	7416-E2	7416-F1	7416-F2
10g	7521-E1P	7521-E2P	7521-F1P	7521-F2P	2kg – 1mg	7415-E1	7415-E2	7415-F1	7415-F2
20g	7520-E1P	7520-E2P	7520-F1P	7520-F2P	2kg – 1g	7414-E1	7414-E2	7414-F1	7414-F2
50g	7518-E1P	7518-E2P	7518-F1P	7518-F2P	5kg – 1mg	7411-E1	7411-E2	7411-F1	7411-F2
100g	7517-E1P	7517-E2P	7517-F1P	7517-F2P	5kg – 1g	7410-E1	7410-E2	7410-F1	7410-F2
200g	7516-E1P	7516-E2P	7516-F1P	7516-F2P	5kg – 1kg	7442-E1	7442-E2	7442-F1	7442-F2
500g	7514-E1P	7514-E2P	7514-F1P	7514-F2P	10kg – 1mg	7409-E1	7409-E2	7409-F1	7409-F2
1kg	7513-E1P	7513-E2P	7513-F1P	7513-F2P	10kg – 1g	7408-E1	7408-E2	7408-F1	7408-F2
2kg	7512-E1P	7512-E2P	7512-F1P	7512-F2P	20kg – 1mg	7407-E1	7407-E2	7407-F1	7407-F2
5kg	7510-E1P	7510-E2P	7510-F1P	7510-F2P	20kg – 1g	7406-E1	7406-E2	7406-F1	7406-F2
10kg	7509-E1P	7509-E2P	7509-F1P	7509-F2P	UKAS / NIST / NVLAP Certificate are available as an option. 选购时可另外附加 UKAS / NIST / NVLAP 校正证书				
20kg	7508-E1P	7508-E2P	7508-F1P	7508-F2P					
50kg	7505-E1	7505-E2	7505-F1	7505-F2P					



香港科学仪器社有限公司
Che Scientific Co. (Hong Kong) Ltd.
Tel: (852) 2481 1323
Fax: (852) 2418 1302

浙江科通仪器有限公司 (绍兴)
Zhejiang Fortune Scientific Co. Ltd. (Shaoxing)
Tel: (86 575) 8812 0560
Fax: (86 575) 8812 0561

浙江科通仪器有限公司 (广州)
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)
Tel: (86 20) 8130 7388
Fax: (86 20) 8108 4925

上海香科仪器贸易有限公司
Shanghai Xiang Ke Scientific Trading Co. Ltd.
Tel: (86 21) 6248 6170
Fax: (86 21) 6248 0647