

WEIFO®

“维和”牌 微电脑防潮箱 “WEIFO” Electronic Auto-Drying Cabinet

台湾制造
Made in Taiwan

Standard Series Electronic Auto-Drying Cabinet 实用型 微电脑防潮箱

数字显示湿度, 20-60%RH 湿度设定, DRY-205 型同时显示温度
另有 Super Dry 型号, 10-60%RH 湿度设定



Economy Series Electronic Auto-Drying Cabinet 经济型 微电脑防潮箱

棒条式显示湿度, 30 / 40 / 50 / 60%RH 湿度设定



MODEL 型号	实用型				经济型	
	DRY-60	DRY-70	DRY-100	DRY-205	KD-70	KD-112
External Dim. (WxDxHmm) 外形尺寸 (阔×深×高 mm)	400x410x560	400x410x660	400x410x890	400x530x1280	400x410x560	400x410x890
Internal Dim. (WxDxHmm) 内部尺寸 (阔×深×高 mm)	398x355x480	398x355x580	398x355x820	398x475x1210	398x355x480	398x355x820
Capacity 有效容积	68L(公升)	82L(公升)	116L(公升)	228L(公升)	70L(公升)	116L(公升)
Door 门	SINGLE (单门)	SINGLE (单门)	DUAL(双门)	DUAL(双门)	SINGLE (单门)	SINGLE (单门)
Shelf 隔板	2 PCS (片)	2 PCS(片)	3 PCS(片)	5 PCS(片)	2 PCS(片)	3 PCS(片)
Weight 重量((约)	15kg (千克)	17kg (千克)	20kg (千克)	30kg (千克)	15kg (千克)	20kg (千克)

Professional Series Electronic Auto-Drying Cabinet 专业型 微电脑防潮箱

数字显示湿度及温 度, 20-60%RH 湿度设定, 另有 Super Dry 型号, 10-60%RH 湿度设定



PROFESSIONAL SERIES 专业型							
Model 型号	DRY-125	DRY-300	DRY-420	DRY-520	DRY-610	DRY-530	DRY-1360
External Dim. (WxDxHmm) 外形尺寸 (阔×深×高 mm)	440x440x930	890x480x930	1180x480x1020	1180x570x1020	1400x600x100 连脚轮	920x660x1190 连脚轮	1200x730x1900 连脚轮
Internal Dim. (WxDxHmm) 内部尺寸 (阔×深×高 mm)	438x385x850	888x425x850	1178x425x900	1178x515x900	1398x545x850	888x610x1070	1198x675x1740
Capacity 有效容积	141L(公升)	320L(公升)	450L(公升)	545L(公升)	646L(公升)	578L(公升)	1402L(公升)
Door 门	SINGLE (单门)	DUAL (双门)	DUAL (双门)	DUAL (双门)	DUAL (双门)	4 PCS DOORS (FRONT&BACK) (前后双门)	6 PCS DOOR (六门)
Shelf 隔板	3 PCS (片)	3 PCS (片)	3 PCS (片)	3 PCS (片)	3 PCS (片)	6 PCS (片)	5 PCS (片)
Weight 重量((约)	25kg (千克)	43kg (千克)	55kg (千克)	60kg (千克)	68kg (千克)	70kg (千克)	120kg (千克)

香港科学仪器社有限公司
Che Scientific Co. (Hong Kong) Ltd.
Tel: (852) 2481 1323
Fax: (852) 2418 1302

浙江科通仪器有限公司(绍兴)
Zhejiang Fortune Scientific Co. Ltd. (Shaoxing)
Tel: (86 575) 8812 0560
Fax: (86 575) 8812 0561

浙江科通仪器有限公司(广州)
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)
Tel: (86 20) 8130 7388
Fax: (86 20) 8108 4925

上海香科仪器贸易有限公司
Shanghai Xiang Ke Scientific Trading Co. Ltd.
Tel: (86 21) 6248 6170
Fax: (86 21) 6248 0647