


Microscope Accessories



MS-1

Model MS-1 Graduated Mechanical Stage. This full size graduated mechanical stage has a range of movement from 26mm to 67mm with a accuracy reading to 0.1mm. the stage is made to fit most microscopes that are predrilled to accept mechanical stages.

HKS-SP-32

Substage Illuminator, Model HKS-SP-32, for use with most standard microscopes that come with mounted in the base or on the arm of the microscope. with 110V/120V AC cord and plug. 110/120V AC - 20W bulb is included. (220V AC also available.)



HKS-SP-33A to HKS-SP-33D

Achromatic Objectives. The objectives below are standard size, coated achromats. They will fit most standard microscopes.

Model

HKS-SP-33A	4X
HKS-SP-33B	10X
HKS-SP-33C	40X(S)
HKS-SP-33D	100X(S), Oil.

Mirrors and Forks. Designed to fit most compact and standard microscopes.

Model


HKS-SP-34A, 40mm diameter (for microscope model HKS--7A)

HKS-SP-34B, 50mm diameter. (for standard microscopes)



HKS-SP-34A to HKS-SP-34B

HKS-SP-3B to 3C



HKS-SP-3

Microscope Bulb. For microscope or microscope series.

Model

HKS-SP-3	Halogen Bulb, 6V 20W
HKS-SP-3B	Tungsten bulb, 110V 20W
HKS-SP-3C	- do - , 220V 20W