

“Globe” Brand, Grade B, Economy Meter, These full size meters offer quality and value. Each meter delivers accuracy of $\pm 2.5\%$. Features such as external zero-adjust calibration and combine with a large 7.5 cm scale. Easy to read inclined meter face, and a fully enclosed plastic case to make them excellent for student use. All meters are in black plastic cases. DC are in moving coil type and DC/AC are in moving iron type.



“Grade B” Voltmeter (Black color case)

Voltmeters (Grade B)

- PA-1 Voltmeter, triple range : 0-3V/15V/300V DC, in 0.1V/0.5V/10V.
- PA-1-1 Voltmeter, triple range : 0-3V/10V/15V. DC, in 0.1V/0.2V/0.5V
- PA-1-1-1 Voltmeter, triple range : 0-3V/15V/30V DC in 0.1V/1V/1V
- PA-2 Voltmeter, range : 0-3V. in 0.1V. DC.
- PA-3 Voltmeter, range : 0-5V. in 0.1V. DC.
- PA-3-1 Voltmeter, range : 0-10V. in 0.2V. DC.
- PA-4 Voltmeter, double range: 0-5V/0-15V. DC. in 0.1V/0.3V.
- PA-4-1 Voltmeter, range : 0-15V DC in 0.3V

Ammeters (Grade B)

- PA-5 Ammeter, triple range : 0-50mA/500mA/5A, in 1mA/10mA/100mA, DC.
- PA-6 Ammeter, double range : 0-1A/0-5A in 20 mA/100mA, DC.
- PA-7 Ammeter, 0-1A. in 20 mA. DC.
- PA-8 Ammeter, 0-3A. in 100 mA. DC
- PA-9 Ammeter, 0-5A. in 100 mA. DC
- PA-10 Ammeter, 0-10A. in 200 mA. DC
- PA-10A Ammeter, but moving iron type, 0-3A in 0.1A. DC/AC.
- PA-10B Ammeter, moving iron type, 0-10A in 200 mA. DC/AC.
- PA-11 Milliammeter, 0-100 mA in 2 mA. DC
- PA-11-1 Milliammeter, 0-1 mA in 0.02mA. DC
- PA-11-2 Milliammeter, 0-200 mA in 4mA. DC.
- PA-11-3 Milliammeter, 0-500 mA. in 10mA, DC.
- PA-11-4 Milliammeter, 0-10 mA. in 0.2mA, DC.
- PA-12 Microammeter, -20 to 100 μ A. in 2 μ A. DC
- PA-13 Microammeter, 0-500 μ A in 10 μ A. DC



“Grade B” Ammeter (Black color case)



“Grade B” Galvanometer (Black color case)

Galvanometers (Grade B)

- PA-14 Sensitive Galvanometer, -35mV/+35mV, in 1mV. DC.
- PA-14-1 Sensitive Galvanometer, -500 μ A/0/+500 μ A \pm 20 μ A
- PA-14-2 Sensitive Galvanometer, -50 μ A/0/+50 μ A. \pm 2 μ A