




HC-1

**Hand Tally Counter, Model HC-1**, this easy-to-use hand counter has four wheels for counting to to 9,999. It is useful for counting anything from pendulum swings to bacteria cultures. Mounted in a chrome-plated metal case.



**Model**

- PA-131** Single Pole, Single Throw
- PA-132** Single Pole, Double Throw,
- PA-133** Double Pole, Single Throw
- PA-134** Double Pole, Double Throw

**Plastic Knife Switch**, these sturdy knife switched feature bras poles and knives, making the switches well-suited in general-purpose classroom and laboratory experiments, all come with binding posts for ease of wire connection.



PA-122-BC and PA-122-BF

**Digital LCD Electronic Clinical Thermometer**, easy-to-read temperature with this inexpensive digital thermometer with LCD. The solid-state electronics are housed in an unbreakable plastic case and batteries are included.

**Model**

- PA-122-BC** readout. in  $^{\circ}\text{J}$ . Resolution: 0.1  $^{\circ}\text{J}$  (range : 32  $^{\circ}\text{J}$  to 42  $^{\circ}\text{J}$ )
- PA-122-BF** readout in  $^{\circ}\text{K}$  Resolution: 0.1  $^{\circ}\text{K}$ . (range : 95  $^{\circ}\text{K}$  to 102  $^{\circ}\text{K}$ )

**Glass Tubing Cutter, Model PA-145**, hardened steel wheel cuts glass tubes up to 38 mm dia. Two V guides align the tubing, wheel is spring loaded to stay open. Zinc alloy.



PA-145