






MICROSCOPE GLASS SLIDES & COVER GLASS

Microscope Glass Slide, made of clear high quality non-corrosive glass, interleaved, ground edges and pre-cleaned. All microscope slides are 25.4 x 76.2 mm. (1" x 3").



Model

- MG-1** Microscope Slides, thickness: 1.0 - 1.2 mm. Box of 72 pcs.
MG-1A - ditto - but 50 x 76 mm. (2" x 3"). Box of 72 pcs.
MG-1D Microscope Slides, thickness: 1.0 - 1.2 mm. Box of 50 pcs.
MG-2 Microscope Slide. but thickness: 0.8-1.0 mm., size: 1" x 3". Box of 72 pcs.
MG-3 Microscope Slides, but frosted in one end on one side. Box of 72 pcs.
MG-4 Microscope Slide, with one cavity. Diam. 15 mm (5/8"), size: 25.4 x 76.2 mm. (1" x 3"). Box of 72 pcs.
MG-4A - ditto - but with two cavity. Box of 72 pcs.
MG-4C - ditto - but thickness 5 mm (1/4"), with diam. 16 mm (3/4"), single 3.25 mm (3/16"). Deep Well. Box of 50 pcs.
MG-4D Microscope Slide. With one cavity diam. 15 mm (5/8")., size: 25.4 x 76.2 mm. (1" x 3"). Box of 50 pcs.
MG-5 Microscope Cover Glass. size: 18x18 mm (3/4" x 3/4"), thickness No. 1 Box of 100 pcs.
MG-6 - ditto - but size: 22 x 22 mm. (1" x 1"). Box of 100 pcs.
MG-7 - ditto - but size: 20 x 20 mm. (7/8" x 7/8"). Box of 100 pcs.

Filter Paper and Lens Tissue Paper

The mostly widely used general purpose paper. Good quality for qualitative analysis. Medium rate. 100 circles per box.



Model

- FP-6** Diam. 7 cm.
FP-1 Diam. 9 cm.
FP-2 Diam. 11 cm.
FP-3 Diam. 12.5 cm.
FP-4 Diam. 15 cm.
FP-5 Diam. 18 cm.

Model

- LP-7** Lens tissue paper. 10 x 15 cm, pack of 100 sheets.
LP-10 Test Paper, litmus, blue, pack of 80 strips.
LP-11 - do -, Red.
LP-12 - do -, Potassium Iodine Starch.
LP-13 - do -, Phenolphthalein.
LP-14 - do -, Ph.