

## Glassware

**Glass Beaker**, low form, with spout, borosilicate glass, with printed graduations for 50 - 1000 ml.




GW-13 to GW-21


<u>Model</u>	<u>Cap. (ml)</u>	<u>O.D. Diam. (mm)</u>	<u>Height (mm)</u>
<b>GW-13</b>	50	46	58
<b>GW-14</b>	100	52	72
<b>GW-15</b>	150	57	80
<b>GW-16</b>	250	69	94
<b>GW-17</b>	400	82	110
<b>GW-18</b>	500	87	115
<b>GW-19</b>	600	95	128
<b>GW-20</b>	1000	110	150
<b>GW-21</b>	2000	133	200

**Erlenmeyer Flask**, conical narrow mouth, borosilicate glass.

<u>Model</u>	<u>Cap.</u>	<u>Neck O.D. Diam.</u>	<u>Body O.D. Diam.</u>	<u>Height</u>
<b>GW-22</b>	250ml	30 mm.	83 mm.	144 mm.
<b>GW-23</b>	500ml	35 mm.	100 mm.	195 mm.
<b>GW-24</b>	1000ml	40 mm.	128 mm.	225 mm.
<b>GW-25</b>	2000ml	53 mm.	170 mm.	285 mm.



GW-22 to GW-25



GW-26 to GW-29

**Boiling Flask**, flat bottom, narrow neck, borosilicate glass.

<u>Model</u>	<u>Cap.</u>	<u>Neck O.D. Diam.</u>	<u>Body O.D. Diam.</u>	<u>Height</u>
<b>GW-26</b>	250ml	25 mm.	88 mm.	160 mm.
<b>GW-27</b>	500ml	30 mm.	110 mm.	200 mm.
<b>GW-28</b>	1000ml	35 mm.	140 mm.	250 mm.
<b>GW-29</b>	2000ml	40 mm.	175 mm.	300 mm.

**Boiling Flask**, round bottom, borosilicate glass.

<u>Model</u>	<u>Cap.</u>	<u>Neck O.D. Diam.</u>	<u>Body O.D.Diam.</u>	<u>Height</u>
GW-30	250ml	25 mm.	88 mm.	165 mm.
GW-31	500ml	30 mm.	110 mm.	210 mm.
GW-32	1000ml	35 mm.	140 mm.	260 mm.
GW-33	2000ml	46 mm.	175 mm.	305 mm.



GW-30 to GW-33



**Volumetric Flask**, soft soda glass, graduated one mark with plastic stopper.

<u>Model</u>	<u>Cap.</u>	<u>Neck O.D. Diam.</u>	<u>Body O.D.Diam.</u>	<u>Height</u>
GW-34	50ml	11 mm.	46 mm.	140 mm.
GW-35	100ml	12 mm.	59 mm.	182 mm.
GW-36	250ml	14 mm.	80 mm.	225 mm.
GW-37	500ml	17 mm.	101 mm.	250 mm.
GW-38	1000ml	19mm.	126 mm.	320 mm.

GW-34 to GW-38

**Graduated Cylinders with spout**, soft soda glass, single, in brown graduation.

<u>Model</u>	<u>Cap.</u>	<u>O.D. Diam.</u>	<u>Height</u>	<u>Graduated interval</u>
GW-5	10 ml	13 mm.	150 mm.	0.2 ml.
GW-6	25 ml	20 mm.	170 mm.	0.5 ml.
GW-7	50 ml	25 mm.	195 mm.	1.0 ml.
GW-8	100 ml	30 mm.	250 mm.	2.0 ml.
GW-9	250 ml	40 mm.	320 mm.	5.0 ml.
GW-10	500 ml	54 mm.	360 mm.	5.0 ml.
GW-11	1000 ml	65 mm.	440 mm.	10.0 ml.
GW-12	2000 ml	90 mm.	500 mm.	20.0 ml



GW-5 to GW-12



**Glass Gas Jar**, soft soda glass with ground on cover plate..

- Model  
**GW-1 50 x 150mm (D. x H.)**  
**GW-2 50 x 200mm (D.x H.)**

GW-1 & GW-2



GW-3

**Glass Alcohol Lamp, Model GW-3**, conical shape, with plastic cap 150 ml but with wick.