

Laboratory Equipment 实验室设备



“李可”牌 离心机
“BOECO” Centrifuge

德国
Germany

SC-8 Compact Centrifuge 小型离心机 (订货编号: BOE 01800-13)

- 标配8位固定角转轴, 最大可放 15ml 玻璃管或锥形管
- 最高转速: 6,000RPM / 3,421RCF
- 转速显示可选 RPM 或 RCF
- 定时器 1sec – 99 mins
- 符合标准: EN 55011 ISM class B, EN 61000-6-1, IEC 61010
- 电源: 220V 50Hz AC

capacity in ml	4	5	7	15	15	blood collection / urine tubes	
dimensions Ø x L mm	10 x 88	12/13 x 75	13 x 100	17 x 100	17 x 120	max 17 x 100	
n=6,000 min ⁻¹							
Rotor included							
Cat.-No.	06305	01054-A	01058	-	-	-	
tubes per rotor	8	8	8	8	4	8	
RCF	2,817	2,697	3,461	3,461		3,461	
radius in mm	70	67	86	86		99	
run-up in sec						26	
run-down in sec. ²⁾						37	



HC-240 Hematocrit Centrifuge 离心机 (订货编号: BOE 01801-13)

- 适用于血球容积比测定(hematocrit determinations)
- 标配转轴可放 24 支毛细管
- 最高转速: 13,000RPM / 16,060RCF
- 转速显示可选 RPM 或 RCF
- 定时器 1sec – 99 mins
- 符合标准: EN / IEC 61326-1 class B / FCC class B, IEC 61010
- 电源: 220V 50Hz AC

capacity	75 µl
Capillaries	1.5 x 75 mm
Rotor and Lid as evaluation disc incl.	
Cat.-No.	-
capillaries per rotor	24
RCF	16,060
radius in mm	85
run-up in sec	8
run-down in sec. ²⁾	10



转轴连盖
(Lid as evaluation disc included)



MC-24 Microlitre Centrifuge 微升离心机 (订货编号: BOE 01203-13)

- 可选不同固定角转轴, 最多可放 24 支 1.5ml / 2.0ml 微量离心管
- 最高转速: 13,300RPM / 17,008RCF
- 转速显示可选 RPM 或 RCF
- 定时器 1sec – 99 mins
- 符合标准: EN / IEC 61326-1 class B / FCC class B, IEC 61010
- 电源: 220V 50Hz AC

capacity in ml	1.5	0.5/0.8	0.2 / 0.4	2.0
dimensions Ø x L mm	11x38	8x30/45	6x18/45	11x38
n=13,300 min ⁻¹				
Cat.-No. 01252	with high speeds we recommend to use adapter No. 02031			
Cat.-No.	-	02023	02024	-
tubes per rotor	12			
RCF	14,041			
radius in mm	71			
run-up in sec.	15			
run-down in sec. ²⁾	16			



Angle rotor,
12-place
Cat. No. 01252
12-place,
n= 13,300 min⁻¹,
max RCF 14,041



capacity in ml	1.5	0.5/0.8	0.2 / 0.4	2.0
dimensions Ø x L mm	11x38	8x30/45	6x18/45	11x38
n=13,300 min ⁻¹				
Cat.-No. 01258	with high speeds we recommend to use adapter No. 02031			
Cat.-No.	-	02023	02024	-
tubes per rotor	18			
RCF	15,030			
radius in mm	76			
run-up in sec.	15			
run-down in sec. ²⁾	16			



Angle rotor,
18-place,
Cat. No. 01258
n= 13,300 min⁻¹,
max RCF 15,030



Angle rotor, 24-place,
Cat. No. 01226
n= 13,300 min⁻¹,
max RCF 17,008

capacity in ml	1.5	0.5/0.8	0.2 / 0.4	2.0
dimensions Ø x L mm	11x38	8x30/45	6x18/45	11x38
n=13,300 min ⁻¹				
Cat.-No. 01226	with high speeds we recommend to use adapter No. 02031			
Cat.-No.	-	02023	02024	-
tubes per rotor	24			
RCF	17,008			
radius in mm	86			
run-up in sec.	15			
run-down in sec. ²⁾	16			

香港科学仪器社有限公司
Che Scientific Co. (Hong Kong) Ltd.
Tel: (852) 2481 1323
Fax: (852) 2418 1302

浙江科通仪器有限公司(绍兴)
Zhejiang Fortune Scientific Co. Ltd. (Shaoxing)
Tel: (86 575) 8812 0560
Fax: (86 575) 8812 0561

浙江科通仪器有限公司(广州)
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)
Tel: (86 20) 8130 7388
Fax: (86 20) 8108 4925

上海香科仪器贸易有限公司
Shanghai Xiang Ke Scientific Trading Co. Ltd.
Tel: (86 21) 6248 6170
Fax: (86 21) 6248 0647