

Culture Tube, with Rubber Lined Screw Cap, Unattached

- Tubes made of KG-33 borosilicate glass
- ASTM E982, Type VI, Class A requirements

订货编号	外径 x 长度(mm)	容量(mL)	包装数量
45066-13100	13 x 100	8	288/CS
45066-16100	16 x 100	12	288/CS
45066-16125	16 x 125	16	288/CS
45066-16150	16 x 150	20	288/CS
45066-20125	20 x 125	25	192/CS
45066-20150	20 x 150	30	192/CS
45066-25150	25 x 150	50	144/CS
45066-25200	25 x 200	70	144/CS
45066-38200	38 x 200	160	24/CS



Culture Tube, Rim top, with Marking Spot

- Tubes made of KG-33 borosilicate glass
- ASTM E982, Type IV requirements

订货编号	外径 x 长度(mm)	容量(mL)	包装数量
45048-1075	10 x 75	4	720/CS
45048-1275	12 x 75	6	720/CS
45048-13100	13 x 100	10	720/CS
45048-15125	15 x 125	16	720/CS
45048-16100	16 x 100	15	576/CS
45048-16125	16 x 125	20	576/CS
45048-16150	16 x 150	24	576/CS
45048-18150	18 x 150	28	576/CS
45048-20150	20 x 150	36	576/CS
45048-25150	25 x 150	55	288/CS
45048-25250	25 x 250	95	144/CS



Culture Tube, with PTFE-Faced Rubber Lined Screw Cap, Unattached

- Tubes made of KG-33 borosilicate glass
- ASTM E982, Type VI, Class B requirements

订货编号	外径 x 长度(mm)	容量(mL)	包装数量
45066A-13100	13 x 100	8	288/CS
45066A-16100	16 x 100	12	288/CS
45066A-16125	16 x 125	16	288/CS
45066A-16150	16 x 150	20	288/CS
45066A-20125	20 x 125	25	192/CS
45066A-20150	20 x 150	30	192/CS
45066A-25150	25 x 150	50	144/CS
45066A-25200	25 x 200	70	144/CS



Culture Tube, Rim top, Disposable

- Tubes made of N-51A borosilicate glass
- ASTM E438, Type I, class B requirements

订货编号	外径 x 长度(mm)	包装数量
73500-650	6 x 50	1000/CS
73500-1075	10 x 75	1000/CS
73500-1275	12 x 75	1000/CS
73500-13100	13 x 100	1000/CS
73500-16100	16 x 100	1000/CS
73500-16125	16 x 125	1000/CS
73500-16150	16 x 150	1000/CS
73500-18150	18 x 150	500/CS
73500-20150	20 x 150	500/CS
73500-25150	25 x 150	500/CS



DUALL® Tissue Grinder, All-Glass

- Combines both conical and cylindrical surfaces to effectively reduce tissue and produce a uniform homogenate.
- This ground glass model is ideal for connective tissue, e.g. muscle, heart, lung
- Designed for a motor drive using a 788000 Stirrer Adapter.

订货编号	管长 x 外径(mm)	容量(mL)	包装数量
885450-0020	80 x 13	1	1/PK
885450-0021	120 x 16	3	1/PK
885450-0022	150 x 18	5	1/PK
885450-0023	175 x 25	15	1/PK
885450-0024	215 x 32	30	1/PK
885450-0025	225 x 38	50	1/PK



Model 885480 with PTFE pestle and glass tube.

KIM-KAP®, Polypropylene Closures for Rim top tubes, Autoclavable, Color Coded

- Other colors are also available: 73662 (Red), 73663 (Green), 73664 (Yellow), 73665 (Blue)

订货编号	适用外径 (mm)	颜色	包装数量
73660-13	13	Natural	1000/CS
73660-16	16	Natural	1000/CS
73660-18	18	Natural	1000/CS
73660-20	20	Natural	500/CS
73660-25	25	Natural	500/CS
73660-38	38	Natural	250/CS



Ampuls

- 12010L - Trimmed Stem
- 12040G - Trimmed Stem, Amber

订货编号	容量(mL)	包装数量
12010L-1	1	1152/CS
12010L-2	2	1152/CS
12010L-5	5	864/CS
12010L-10	10	576/CS
12010L-20	20	288/CS
12040G-1	1	144/CS
12040G-2	2	144/CS



12010L
标准



12040G
茶色

Potter-Elvehjem Tissue Grinder, All-Glass

- Designed for a motor drive using a 788000 Stirrer Adapter.

订货编号	管长 x 外径(mm)	容量(mL)	包装数量
885500-0019	145 x 5	1	1/PK
885500-0021	205 x 6	5	1/PK
885500-0022	215 x 8	8	1/PK
885500-0023	265 x 10	17	1/PK
885500-0024	315 x 10	45	1/PK



Model 886000 with PTFE pestle and glass tube.

Laboratory Equipment 实验室设备



“金布尔”牌 实验室玻璃设备
“Kimble” Laboratory Glassware

美国
U.S.A.

ULTRA-WARE® Microfiltration Assembly, 47mm, with Support

Fritted Glass Support:

- Used for general filtration.

订货编号	漏斗容量 (mL)	过滤器直径 (mm)
953755-0000	300	47

Component Parts:

订货编号	描述
953751-0000	1. Glass Funnel, 47mm, 300mL
953753-0000	2. Anodized Aluminum Clamp, 47mm
953752-0001	3. Fritted Glass Support, Base, 47mm
953763-0801	4. #8 Silicone Stopper, 9/16" (14.3mm) Hole




Stainless Steel Support:

- For filtering viscous or proteinaceous solutions to give the maximum flow rate.

订货编号	漏斗容量 (mL)	过滤器直径 (mm)
953805-0000	300	47

Component Parts:

订货编号	描述
953751-0000	1. Glass Funnel, 47mm, 300mL
953753-0000	2. Anodized Aluminum Clamp, 47mm
953810-0000	3. Stainless Steel Support Screen
953811-0000	4. PTFE Support Screen Gasket, Pkg./25, (1 Supplied)
953801-0001	5. Glass Support Base, 47mm
953763-0801	6. #8 Silicone Stopper, 9/16" Hole



Microfiltration Accessories

订货编号	描述
953761-0000	Glass Funnel, 47mm, 100mL Capacity
953771-0000	Glass Funnel, 47mm, 500mL Capacity
953781-0000	Glass Funnel, 47mm, 1000mL Capacity
953760-0000	Filtration Flask, 1L, #8 Stopper Joint, 1/4" Hose
953760-1002	Filtration Flask, Safety Coated, 1L, #8 Stopper Joint, 1/4" Hose



ROCKER Filtration Manifold

(Made in Taiwan)

- Each branch is individually controlled

订货编号	描述
167200-43	3-position, aluminum
167110-23	3-position, SUS316
167110-26	6-position, SUS316



Volumetric Flask, Class A, w/ PTFE Stopper

- With a graduation ring blasted on the neck.
- The 5 mL and larger sizes are designed from ASTM Specification E288, Class A unserialized requirements.
- Replacement stopper is 41901R.

订货编号	容量 (mL)	公差 (±mL)	塞子大小	包装数量
28014F-2	2	0.015	8	12/CS
28014F-5	5	0.02	9	12/CS
28014F-10	10	0.02	9	12/CS
28014F-25	25	0.03	9	12/CS
28014F-50	50	0.05	9	12/CS
28014F-100	100	0.08	13	12/CS
28014F-200	200	0.1	16	12/CS
28014F-250	250	0.12	16	12/CS
28014F-500	500	0.2	19	12/CS
28014F-1000	1000	0.3	22	6/CS
28014F-2000	2000	0.5	27	4/CS



28014F-100

Separatory Funnel, Squibb, Pear-shaped, w/ PTFE Stopcock, Plastic Stopper

- Designed from ASTM Specification E1096, type 4 requirements.
- Replacement stopcock plug is a 41500F.
- Replacement stopper is 28160R.

订货编号	容量 (mL)	活塞大小	塞子大小	包装数量
29049F-125	125	2	22	4/CS
29049F-250	250	4	22	4/CS
29049F-500	500	4	27	4/CS
29049F-1000	500	4	27	2/CS
29049F-2000	1000	6	38	2/CS



29049F-125

B.O.D. Bottles with Glass Stopper

- 300mL Capacity
- For Method 5210-B: 5-day B.O.D. Test
- Conform to ASTM E438, Type I, Class B


订货编号	数列	包装数量
15070-00	Un-numbered	24/CS
15070-01	1-24	
15070-25	25-48	

Replacement Bottle:


15070L-99	Un-numbered	24/CS
-----------	-------------	-------

Snap-Lid holds stoppers secure over the flared mouth

订货编号	包装数量
15070L-99	Snap-Lid 50/CS



Snap-Lid Securely hold stoppers



Laboratory Equipment 实验室设备



“SciLabware” 实验室玻璃设备
“SciLabware” Laboratory Glassware

英国
England

PYREX Burettes, PTFE key, class AS

订货编号	容量 (mL)	刻度 (mL)	误差 (±mL)
3295/12	10	0.02	0.02
3295/14	25	0.10	0.05
3295/16	50	0.10	0.05
3295/18	100	0.20	0.10

PYREX Burettes, PTFE key, class B

订货编号	容量 (mL)	刻度 (mL)	误差 (±mL)
3295/02	10	0.02	0.05
3295/04	25	0.10	0.10
3295/06	50	0.10	0.10
3295/08	100	0.20	0.20



3295

PYREX Cylinders, with spout, class A

订货编号	容量 (mL)	刻度 (mL)	误差 (±mL)	包装数量
3201/02	5	0.1	0.05	2
3201/04	10	0.2	0.10	2
3201/06	25	0.5	0.25	2
3201/08	50	1.0	0.50	2
3201/10	100	1.0	0.50	2
3201/12	250	2.0	1.00	2
3201/14	500	5.0	2.00	2
3201/16	1000	10.0	5.0	2
3201/18	2000	20.0	10.00	1



PYREX Cylinders, white graduations, class B

订货编号	容量 (mL)	刻度 (mL)	误差 (±mL)	包装数量
3200/02MP	5	0.1	0.10	2
3200/04MP	10	0.2	0.20	2
3200/06MP	25	0.5	0.50	2
3200/08MP	50	1.0	1.00	2
3200/10MP	100	1.0	1.00	2
3200/12MP	250	2.0	2.00	2
3200/14MP	500	5.0	5.00	2
3200/16MP	1000	10.0	10.00	2
3200/18MP	2000	20.0	20.00	1



MBL Burettes, economy, PTFE key, class B

➢ BRB074P, BRB077P with fluted funnel tops

订货编号	容量 (mL)	刻度 (mL)	误差 (±mL)
BRB070P	10	0.02	0.05
BRB074P	25	0.10	0.10
BRB075P	25	0.10	0.10
BRB076P	50	0.10	0.10
BRB077P	50	0.10	0.10



Funnel tops

BRB070P

PYREX Crucibles, Gooch

订货编号	容量 (mL)	烧结碟直径 (mm)	孔隙率	包装数量
3650/02M	10	20	1	5
3650/04M	10	20	2	5
3650/06M	10	20	3	5
3650/08M	10	20	4	5
3660/02M	30	30	1	5
3660/04M	30	30	2	5
3660/06M	30	30	3	5
3660/08M	30	30	4	5
3660/10M	30	30	5	5
3680/02M	60	40	1	5
3680/04M	60	40	2	5
3680/06M	60	40	3	5
3680/08M	60	40	4	5
3680/10M	60	40	5	5



- Ideal for drying precipitates to constant mass at 110°C
- Also suitable for ignition work up to 450°C

PYREX Evaporating Dishes

订货编号	容量 (mL)	直径 (mm)	包装数量
1450/01D	15	50	10
1450/03D	45	60	10
1450/04D	90	80	10
1450/06D	170	95	10
1450/10D	320	115	10
1450/14D	600	140	10



QUICKFIT Flasks, round bottom, short neck

订货编号	容量 (mL)	卡口 (mm)	包装数量
FR250/1S	250	14/23	2
FR250/2S	250	19/26	5
FR250/3S	250	24/29	5
FR250/4S	250	29/32	5
FR250/5S	250	34/35	5
FR500/2S	500	19/26	5
FR500/3S	500	24/29	5
FR500/4S	500	29/32	5
FR500/5S	500	34/35	5
FR1L/2S	1000	19/26	2
FR1L/3S	1000	24/29	2
FR1L/4S	1000	29/32	2
FR1L/5S	1000	34/35	2
FR1L/7S	1000	45/40	1
FR2L/3US	2000	24/29	1
FR2L/3S	2000	24/29	1
FR2L/4S	2000	29/32	1
FR2S/5S	2000	34/35	1
FR3L/3S	3000	24/29	1
FR3L/5S	3000	34/35	1
FR4L/7S	4000	45/40	1
FR5L/4S	5000	29/32	1
FR5L/5US	5000	34/35	1
FR5L/9S	5000	55/44	1
FR10L/5S	10000	34/35	1

- With highest quality Quickfit® ground glass socket
- Manufacturer from borosilicate glass for corrosion-free performance and resistance to thermal shock
- Body diameter and wall thickness comply with relevant ISO 1773 and ISO 4797 requirements



香港科学仪器社有限公司
Che Scientific Co. (Hong Kong) Ltd.
Tel: (852) 2481 1323
Fax: (852) 2418 1302

浙江科通仪器有限公司(绍兴)
Zhejiang Fortune Scientific Co. Ltd. (Shaoxing)
Tel: (86 575) 8812 0560
Fax: (86 575) 8812 0561

浙江科通仪器有限公司(广州)
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)
Tel: (86 20) 8130 7388
Fax: (86 20) 8108 4925

上海香科仪器贸易有限公司
Shanghai Xiang Ke Scientific Trading Co. Ltd.
Tel: (86 21) 6248 6170
Fax: (86 21) 6248 0647


Laboratory Equipment 实验室设备



“亨可”牌 实验室玻璃设备
“BOECO” Laboratory Glassware


德国
Germany

Beakers, low form / tall form, with spout, graduated, boro 3.3




容量(mL)	矮杯 Low Form	高杯 Tall Form	包装数量
25	BOE 5010614	---	10
50	BOE 5010617	BOE 5011617	10
100	BOE 5010624	BOE 5011624	10
150	BOE 5010629	BOE 5011629	10
250	BOE 5010636	BOE 5011636	10
400	BOE 5010641	BOE 5011641	10
600	BOE 5010648	BOE 5011648	10
800	BOE 5010653	BOE 5011653	10
1000	BOE 5010654	BOE 5011654	10
2000	BOE 5010663	---	4
3000	BOE 5010668	---	4
5000	BOE 5010673	---	2

Conical Flasks, narrow neck / wide neck, flange-edged, graduated, boro 3.3



容量(mL)	窄口 Narrow Neck	阔口 Wide Neck	包装数量
25	BOE 5021614	BOE 5022614	10
50	BOE 5021617	BOE 5022617	10
100	BOE 5021624	BOE 5022624	10
200	BOE 5021632	BOE 5022632	10
250	BOE 5021636	BOE 5022636	10
300	BOE 5021639	BOE 5022639	10
500	BOE 5021644	BOE 5022644	10
1000	BOE 5021654	BOE 5022654	10
2000	BOE 5021663	BOE 5022663	6
3000	BOE 5021668	---	6
5000	BOE 5021673	---	2


Round / Flat Bottom Flasks, narrow neck, with beaded rim, boro 3.3



容量(mL)	圆底 Round Bottom	平底 Flat Bottom	包装数量
250	BOE 5072136	BOE 5071136	10
500	BOE 5072144	BOE 5071144	10

Vacuum Desiccators / Desiccators, with plastic knob, boro 3.3


- Without porcelain Plate
- Vacuum desiccators with stopcock



直径(mm)	真空活塞 Vacuum Model	普通型号 General Model	包装数量
100	GLA 415262112	GLA 415262111	1
150	GLA 415262162	GLA 415262161	1
200	GLA 415262212	GLA 415262211	1
250	GLA 415262262	GLA 415262261	1
300	GLA 415262303	GLA 415262302	1

Desiccators Plate, porcelain


直径(mm)	型号	包装数量
90	GLA 9180000358	1
140	GLA 9180000354	1
190	GLA 9180000355	1
240	GLA 9180000356	1
280	GLA 9180000357	1



Filtering Flasks, with glass hose connection, boro 3.3

- For vacuum use

容量(mL)	订货编号	包装数量
100	BOE 5120124	10
250	BOE 5120136	10
500	BOE 5120144	8
1000	BOE 5120154	8
2000	BOE 5120163	6



Laboratory Bottles

“BOECO GERMANY” boro 3.3., graduated.

- Complete with blue crew cap and pouring ring of PP.
- Resistant to steam sterilization up to 140°C

容量(mL)	订货编号	包装数量
100	BOE5080245	10
250	BOE5080365	10
500	BOE5080445	10
1000	BOE5080545	10
2000	BOE5080635	10
5000	BOE5080735	6
10000	BOE5080865	1
20000	BOE5080915	1



Laboratory Equipment 实验室设备

实验室玻璃设备, 塑料用品 Laboratory Glassware, Plasticware

AZLON Beakers, straight sided, printed graduations, polypropylene

订货编号	容量 (mL)	刻度 (mL)	包装数量
BDA208P	10	1.0	10
BDA212P	25	2.0	10
BDA218P	50	5.0	10
BDA224P	100	10.0	10
BDA230P	250	25.0	10
BDA234P	400	50.0	10
BDA236P	500	50.0	10
BDA238P	600	50.0	5
BDA240P	1000	100.0	5
BDA242P	2000	200.0	1



AZLON Bottles, round, wide neck, polypropylene

订货编号	容量(mL)	包装数量
BWP0030PN	30	10
BWP0060PN	60	10
BLP0100P	100	10
BWP0100PN	100	10
BLP0125P	125	10
BWP0125PN	125	10
BLP0150P	150	10
BWP0150PN	150	10
BWP0250PN	250	10
BWP0500PN	500	10
BWP1000PN	1000	5
BWP2000PN	2000	5



AZLON Wash Bottles, wide neck, LDPE

订货编号	容量 (mL)	瓶盖颜色	包装数量
WGW230P	250	白	5
WGW231P	250	红	5
WGW232P	250	黄	5
WGW233P	250	绿	5
WGW234P	250	蓝	5
WGW236P	500	白	5
WGW237P	500	红	5
WGW238P	500	黄	5
WGW239P	500	绿	5
WGW235P	500	蓝	5
WGW443P	1000	白	5
WGW445P	1000	红	5
WGW444P	1000	黄	5
WGW447P	1000	绿	5
WGW446P	1000	蓝	5



AZLON Bottles, round, wide neck, amber, polypropylene

订货编号	容量 (mL)	包装数量
BWP0060AP	60	10
BWP0150AP	150	10
BWP0250AP	250	10
BWP0500AP	500	10
BWP1000AP	1000	5
BWP2000AP	2000	5



VITRI Graduated Pitcher

订货编号	容量 (mL)	刻度 (mL)	包装数量
441-081	500	10	12
442-081	1000	10	6
443-081	2000	20	6
444-081	3000	50	6
445-081	5000	100	6



AZLON Cylinders, moulded graduations, class B, polypropylene

订货编号	容量 (mL)	刻度 (mL)	误差 (±mL)	包装数量
CP0010P	10	0.2	0.2	10
CP0025P	25	0.5	0.5	10
CP0050P	50	1	1.00	10
CP0100P	100	1	1.00	5
CP0250P	250	2	2.00	1
CP0500P	500	5	5.00	1
CP1000P	1000	10	10.00	1
CP2000P	2000	20	20.00	1



LMS Volumetric Pipettes, class B

订货编号	容量(mL)	误差(+/-mL)	包装数量
63331401	0.5	0.010	50
63331402	1	0.015	50
63331403	2	0.020	40
63331404	3	0.030	30
63331405	4	0.030	10
63331406	5	0.030	10
63331407	6	0.040	10
63331408	10	0.040	10
63331409	15	0.060	25
63331410	20	0.060	25
63331411	25	0.060	25
63331412	30	0.100	10
63331413	40	0.100	10
63331414	50	0.100	10
63331415	100	0.160	6



LMS Graduated Pipettes, class B

订货编号	容量 (mL)	刻度 (mL)	包装数量
63336151	0.1	0.001	100
63336152	0.2	0.002	100
63336153	0.5	0.01	100
63336154	1	0.01	100
63336155	2	0.02	100
63336156	5	0.05	100
63336157	5	0.1	100
63336158	10	0.1	100
63336159	20	0.1	100
63336160	25	0.1	100
63336181	50	0.2	100

