

KDS100
Single-Syringe Infusion Pump
单针筒注射泵



KDS200
Two-Syringe Infusion Pump
双针筒注射泵



KDS230
Multi-Syringe Infusion /
Withdrawal Pump
多针筒注射 / 提取泵



- 高分辨率触摸式显示屏
- 实时显示
- 可预设程序
- 可同时进行多人多工作操作
- 有防溢槽
- 通过国际认证 UL, CSA, ETL, CE, CB Scheme, RoHs, WEEE

Legato 200
Multi-Syringe Infusion /
Withdrawal Pump
多针筒注射 / 提取泵



Pump Specifications (All pumps available in CE mark version)

Model 型号	KDS100	KDS200	KDS230	Legato 200
Programmable Model 程式型号	---	KDS200P	KDS230P	---
Syringes Max. No. and Size 针筒 最多数量及尺寸	One 10 µl – 60ml each	Two 10 µl – 140ml each Six 20ml – 60ml each Four 100ml – 140ml each	Ten 10 µl – 10ml each	Two 0.5 µl – 140ml each
Min. Flow Rate 最小流率	1 µl/hr (10 µl syringe)	0.01 µl/hr (10 µl syringe)	0.01 µl/hr (10 µl syringe)	0.3 µl/hr (0.5 µl syringe)
Max. Flow Rate 最大流率				
(10ml syringe)	127 ml/hr	21 ml/min	21 ml/min	---
(60ml syringe)	519 ml/hr	86 ml/min	---	---
(140ml syringe)	---	145 ml/min	---	220.97 ml/min
Linear Force 线性压力	20 lb (9 kg)	40 lb (18 kg)	40 lb (18 kg)	75 lb (34 kg)
Accuracy 精确度	+/- < 1%	+/- < 1%	+/- < 1%	+/- < 0.35%
Weight 重量	4.5 lb (2 kg)	9 lb (4 kg)	9.5 lb (4.3 kg)	4.9 Kg (10.75 lb)