

Laboratory Equipment 实验室设备



“Clifton”牌 電熱恆溫水浴箱
“Clifton” Water Baths

英国
England

NE1 系列 恒温水浴箱 Unstirred Thermostatic Water Baths

- 使用旋钮控制温度，每格为 5°C
- 附超温保护

Clifton



NE1-14



NE1-22



NE1-56

型号	NE1-4	NE1-8	NE1-14	NE1-22	NE1-28	NE1-56
容量	4L	8L	14L	22L	28L	56L
本体尺寸 (W x D x H) mm	332x185x290	332x270x290	361x332x290	538x332x290	538x332x290	660x540x290
浴箱尺寸 (W x D x H) mm	300x150x150	300x240x150	325x300x150	500x300x150	500x300x200	620x500x200
温度控制范围	室温 +5°C 至 100°C (不适用于连续沸腾)					
发热线瓦特	400W	800W	1000W	1500W	1500W	2500W
电压	230V 50Hz AC					

NE2D 系列 数显恒温水浴箱 Unstirred Digital Water Baths

- PID 控制器
- 数字显示(LED)，最少读数 0.1°C
- 附超温警报(设定值+/-4°C)及超温保护



NE2-8D



NE2-28D



NE2-56D

型号	NE2-4D	NE2-8D	NE2-14D	NE2-22D	NE2-28D	NE2-56D
容量	4L	8L	14L	22L	28L	56L
本体尺寸 (W x D x H) mm	332x185x290	332x270x290	361x332x290	538x332x290	538x332x290	660x540x290
浴箱尺寸 (W x D x H) mm	300x150x150	300x240x150	325x300x150	500x300x150	500x300x200	620x500x200
温度控制范围	室温+5°C 至 99°C					
发热线瓦特	400W	800W	1000W	1500W	1500W	2500W
电压	230V 50Hz AC					

配件 Accessories

- 可另购不同盖子



不锈钢三角斜盖



不锈钢同心环盖

适用于水浴大小	4L	8L	14L	22L	28L	56L
不锈钢三角斜盖	SL1-4	SL1-8	SL1-14	SL1-22	SL1-22	SL1-56
不锈钢同心环盖	SLR1-4 (2 x 105mm)	SLR1-8 (4 x 83mm)	SLR1-14 (4 x 105mm)	SLR1-22 (6 x 105mm)	SLR1-22 (6 x 105mm)	-----

香港科学仪器社有限公司
Che Scientific Co. (Hong Kong) Ltd.
Tel: (852) 2481 1323
Fax: (852) 2418 1302

浙江科通仪器有限公司(绍兴)
Zhejiang Fortune Scientific Co. Ltd. (Shaoxing)
Tel: (86 575) 8812 0560
Fax: (86 575) 8812 0561

浙江科通仪器有限公司(广州)
Zhejiang Fortune Scientific Co. Ltd. (Guangzhou)
Tel: (86 20) 8130 7388
Fax: (86 20) 8108 4925

上海香科仪器贸易有限公司
Shanghai Xiang Ke Scientific Trading Co. Ltd.
Tel: (86 21) 6248 6170
Fax: (86 21) 6248 0647